
Recombinant mouse ANGPTL8 protein, His

Catalog Number: bs-49019P

Concentration: >0.5 mg/ml

AA Seq: 16-198/198

Predicted MW: 20

Tags: His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 10mM Tris-HCl with 8M Urea.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: ANGPTL8 (Angiopoietin Like 8) is a Protein Coding gene. Diseases associated with ANGPTL8 include Breast Angiosarcoma and Breast Sarcoma. Among its related pathways are Lipoprotein metabolism and Metabolism.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.