

Recombinant human PGA4 protein, C-His (HEK293)

Catalog Number: bs-43144P

Concentration: >0.5mg/ml

AA Seq: 1-388/388

Predicted MW: 41.7

Tags: C-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >0.5% as determined by SDS-PAGE

Purification: AC


Form: Lyophilized or Liquid

Storage: PBS (pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Pepsin is one of the main proteolytic enzymes secreted by the gastric mucosa. Pepsin consists of a single polypeptide chain and arises from its precursor, pepsinogen, by removal of a 41 amino acid segment from the N-terminus. Pepsinogen is synthesized in the stomach lining, and hydrochloric acid, also produced by the gastric mucosa, is necessary to convert the inactive enzyme and to maintain the optimum acidity (pH 1-3) for pepsin function. Pepsin is particularly effective in cleaving peptide bonds involving aromatic amino acids. Pepsin shows extremely broad specificity, and although bonds involving phenylalanine and leucine are preferred, many others are also cleaved to some extent. The amino acid composition of Pepsin C differs from those of pepsinogen and pepsin especially in the content of basic amino acids, glutamic acid, aspartic acid, leucine and isoleucine.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.