
Recombinant mouse H2-Q10 protein

Catalog Number: bs-49099P

Concentration: >0.5 mg/ml

AA Seq: 25-310/325

Predicted MW: 33

Tags: Tag free

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.