

---

## **Recombinant Candidatus Liberibacter solanacearum protein, Trx & S & His**

Catalog Number: bs-49103P

Concentration: >0.5 mg/ml

Species: Candidatus Liberibacter solanacearum

AA Seq: 23-184/184

Predicted MW: 36.3

Tags: N-Trx & S & His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.