

---

## Recombinant HBeAg protein

Catalog Number: bs-1662P

Concentration: >0.5 mg/ml

AA Seq: 1-212/212

Predicted MW: 21

Detected MW: 25 kDa

Tags: Tag free

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 2M Urea with 50mM PBS (pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background:

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.