
DSS1 Antibody Blocking Peptide

Catalog Number: bs-14432P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: SEM1 (SEM1, 26S Proteasome Complex Subunit) is a Protein Coding gene. Diseases associated with SEM1 include Split Foot and Split Hand-Split Foot Malformation. Among its related pathways are RET signaling and Transport of glucose and other sugars, bile salts and organic acids, metal ions and amine compounds.