
ZNF32 Rabbit pAb

Catalog Number: bs-7219R

Target Protein: ZNF32

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog, Horse, Rhesus monkey, Gorilla, Chinese Hamster, Orangutan)

Predicted MW: 31 kDa

Entrez Gene: 7580

Swiss Prot: P17041

Source: KLH conjugated synthetic peptide derived from human ZNF32: 121-220/273.

Purification: affinity purified by Protein A

Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

PRODUCT SPECIFIC PUBLICATIONS

[IF=5.834] Qin Li. et al. ZNF32 prevents the activation of cancer-associated fibroblasts through negative regulation of TGFB1 transcription in breast cancer. FASEB J. 2023 Mar;37(4):e22837 IHC ; Human . 36934389