

bs-7219R

[Primary Antibody]

ZNF32 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000) Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog, Horse, Rhesus monkey, Gorilla, Chinese Hamster, Orangutan) Predicted MW.: 31 kDa Subcellular Location: Nucleus
Clonality: Polyclonal		
GeneID: 7580	SWISS: P17041	
Target: ZNF32		
Immunogen: KLH conjugated synthetic peptide derived from human ZNF32: 121-220/273.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: Preservative: 0.02% Proclin300, Constituents: 1% BSA, 0.01M PBS, pH7.4. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=5.834]** Qin Li. et al. ZNF32 prevents the activation of cancer-associated fibroblasts through negative regulation of TGFB1 transcription in breast cancer. FASEB J. 2023 Mar;37(4):e22837 IHC ;Human. 36934389