

www.bioss.com.cn sales@bioss.com.cn techsupport@bioss.com.cn 400-901-9800

Recombinant human IFN Gamma protein, N-His Tag

Catalog Number: bs-41493P

AA Seq: 24-161/161

Predicted MW: 17

Tags: C-His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl, 150mM NaCl (pH8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Mammalian Interferon gamma is mainly produced by Tlymphocytes and NK cells. It is a

pleiotropic cytokine involved in the regulation of nearly all phases of immune and inflammatory responses, including the activation, growth and differentiation of T cell, B cells, macrophages, NK cells and other cell types such as endothelial cells and fibroblasts. It has weak antiviral and antiproliferative activity, and poteniates the antiviral and anti tumor

effects of IFN alpha / beta (type I interferon). It is upregulated by IL2, FGF basic, EGF and

downregulated by vitamin D3 or DMN. Labile at pH 2.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.