
BLK Rabbit pAb

Catalog Number: bs-6646R

Target Protein: BLK

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: WB (1:500-2000), IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat)

Predicted MW: 58 kDa

Entrez Gene: 640

Swiss Prot: P51451

Source: KLH conjugated synthetic peptide derived from human BLK/MODY11: 231-330/505.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Modulator of beta-cells function, acting through the up-regulation of PDX1 and NKX6-1 and consequent stimulation of insulin secretion in response to glucose.