
Recombinant Bovine Alpha s2 Casein, N-His Tag

Catalog Number: bs-42026P

Concentration: >0.5 mg/ml

AA Seq: 16-222/222

Predicted MW: 30

Tags: N-His Tag

Activity: Not tested

Endotoxin: Not analyzed

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: PBS (pH7.4) with 500mM NaCl.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Important role in the capacity of milk to transport calcium phosphate.

Casocidin-I inhibits the growth of E.coli and S.carnosus.