
Recombinant IcHV-1 ORF53 protein

Catalog Number: bs-49028P

Concentration: >0.5 mg/ml

AA Seq: 1-308/308

Predicted MW: 35

Tags: Tag free

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 10mM Tris-HCl with 8M Urea and 10% Glycerol.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.