

**bs-1530R**

**[ Primary Antibody ]**

## Slc4a7 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

|   |  |
|---|--|
| <p><b>Host:</b> Rabbit</p> <p><b>Clonality:</b> Polyclonal</p> <p><b>Target:</b> Slc4a7</p> <p><b>Immunogen:</b> KLH conjugated synthetic peptide derived from rat Slc4a7 isoform 2: 1-100/1254.</p> <p><b>Purification:</b> affinity purified by Protein A</p> <p><b>Concentration:</b> 1mg/ml</p> <p><b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p><b>Background:</b> Electroneutral sodium- and bicarbonate-dependent cotransporter with a Na(+):HCO<sub>3</sub>(-) 1:1 stoichiometry. Regulates intracellular pH and may play a role in bicarbonate salvage in secretory epithelia. May also have an associated sodium channel activity.</p> | <p><b>Isotype:</b> IgG</p> <p><b>Applications:</b> <b>ELISA</b> (1:5000-10000)</p> <p><b>Reactivity:</b> (predicted: Human, Mouse, Rat)</p> <p><b>Predicted MW.:</b> 140 kDa</p> <p><b>Subcellular Location:</b> Cell membrane</p> |
|---|--|