
Recombinant JSRV-ENV protein, His

Catalog Number: bs-41046P

Concentration: >0.5 mg/ml

AA Seq: 85-378/615

Predicted MW: 33

Detected MW: 38 kDa

Tags: His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0) with 8M Urea.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: bs-0241P is one synthetic peptide derived from HIV1 Surface protein gp120.

The attachment of enveloped viruses to cells and the fusion of viral and cellular membranes are critical early events in the HIV viral infection. This process is mediated by envelope glycoproteins (gp) on the surface of the virus. The human immunodeficiency virus type 1 (HIV-1) envelope glycoprotein, gp160, is proteolytically cleaved into gp120 and gp41, which remain noncovalently associated with one another. gp120 is one of the proteins that forms the envelope of HIV. gp120 projects from the surface of HIV and binds to the CD4 molecule on helper T cells. gp120 has been a logical experimental HIV vaccine because the outer envelope is the first part of the virus that encounters antibody. gp41 is embedded in the outer envelope of HIV that anchors gp120. gp41 also plays a key role in HIV's infection of CD4+ T cells by facilitating the fusion of the viral and cell membranes. The nomenclature of the gp proteins describes their respective molecular masses (e.g., gp160, gp120, gp41).

VALIDATION IMAGES



The purity of the protein is greater than 90% as

determined by reducing SDS-PAGE.