

bs-1899R

[Primary Antibody]

SPATA17 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 128153</p> <p>Target: SPATA17</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human SRG11: 101-200/361.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Strongly expressed in adult testis but weakly expressed in the spleen and thymus. Strongly expressed in round and elongating spermatids, and weakly or not expressed in spermatozoa.</p>	<p>Isotype: IgG</p> <p>SWISS: Q96L03</p> <p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Cow, Dog, Horse)</p> <p>Predicted MW.: 43 kDa</p> <p>Subcellular Location: Nucleus</p>
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------