

---

## Recombinant MeV F protein, His

Catalog Number: bs-42168P

Concentration: >1mg/ml

AA Seq: 24-494/550

Predicted MW: 64.9

Tags: N-His

Activity: Not tested

Purity: >90% as determined by SDS-PAGE

Form: Lyophilized or Liquid

Storage: 20mM Tris (pH8.0) with 8M Urea 0.5M L-Arg.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: No data available.

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.