

**bs-16037R**

**[ Primary Antibody ]**

## **FBX47 Rabbit pAb**



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)  <b>Reactivity:</b> (predicted: Human, Mouse, Rat, Sheep, Cow, Dog)  <b>Predicted MW.:</b> 52 kDa
<b>Clonality:</b> Polyclonal		
<b>GeneID:</b> 494188	<b>SWISS:</b> Q5MNV8	
<b>Target:</b> FBX47		
<b>Immunogen:</b> KLH conjugated synthetic peptide derived from human FBX47: 21-120/452.		
<b>Purification:</b> affinity purified by Protein A		
<b>Concentration:</b> 1mg/ml		
<b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
<b>Background:</b> Probably recognizes and binds to some phosphorylated proteins and promotes their ubiquitination and degradation.		