

bs-16040R

[Primary Antibody]

FBXL22 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 283807</p> <p>Target: FBXL22</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human FBXL22: 161-247/247.</p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: This gene encodes a member of the F-box protein family. This F-box protein interacts with S-phase kinase-associated protein 1A and cullin in order to form SCF complexes which function as ubiquitin ligases.[provided by RefSeq, Sep 2010]</p>	<p>Isotype: IgG</p> <p>SWISS: Q6P050</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500)</p> <p>Reactivity: (predicted: Human)</p> <p>Predicted MW.: 27 kDa</p> <p>Subcellular Location: Cytoplasm ,Nucleus</p>
--	--	---