
NTF4 Antibody Blocking Peptide

Catalog Number: bs-0158P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Neurotrophin 4 (NT4; synonymous with NT5) belongs to the NGF beta family and is a survival factor for peripheral sensory sympathetic neurons. The expression of this gene is ubiquitous and less influenced by environmental signals. NT4 levels are highest in the prostate, with lower levels in thymus, placenta, and skeletal muscle. NT4 is also expressed in embryonic and adult tissues.