
T3-KLH

Catalog Number: bs-0339PK

AA Seq: Coupling protein

Activity: No

Form: Liquid

Storage: 0.01M PBS(pH7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: T3(Triiodothyroxine)although only a minority of the total thyroxine, but the T3 is very important to check thyroid function. In the biochemical function of T4 than strong, and in the case of hyperthyroidism, T3 concentration in the blood are often marked increase (or even up to 5 times the normal).

CAS Number: 51-48-9

Empirical Formula (Hill Notation): C₁₅H₁₂I₃NO₄

Molecular Weight: 650.97

Beilstein Registry Number: 2710227

EC Number: 229-999-3

MDL number: MFCD00002593