
Orai1 Rabbit pAb

Catalog Number: bs-18019R

Target Protein: Orai1

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Polyclonal

Isotype: IgG

Applications: IHC-P (1:100-500), IHC-F (1:100-500), IF (1:100-500), ICC/IF (1:100-500), ELISA (1:5000-10000)

Reactivity: (predicted:Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog, Horse)

Predicted MW: 33 kDa

Entrez Gene: 84876

Swiss Prot: Q96D31

Source: KLH conjugated synthetic peptide derived from human Orai1: 151-250/301.

Purification: affinity purified by Protein A

Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Calcium release activated calcium channel protein 1; Calcium release activated calcium modulator 1; Calcium release-activated calcium channel protein 1; CRACM 1; CRACM1; CRACM1; CRACM1_HUMAN; FLJ14466; Orai 1; ORAI calcium release activated calcium modulator 1; orai1; ORAT1; Protein orai 1; Protein orai-1; TMEM 142A; TMEM142A; Transmembrane protein 142A.