

bs-18019R

[Primary Antibody]

Orai1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

<p>Host: Rabbit</p> <p>Clonality: Polyclonal</p> <p>GeneID: 84876</p> <p>Target: Orai1</p> <p>Immunogen: KLH conjugated synthetic peptide derived from human Orai1: 151-250/301. < Extracellular ></p> <p>Purification: affinity purified by Protein A</p> <p>Concentration: 1mg/ml</p> <p>Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.</p> <p>Background: Calcium release activated calcium channel protein 1; Calcium release activated calcium modulator 1; Calcium release-activated calcium channel protein 1; CRACM 1; CRACM1; CRCM1; CRCM1_HUMAN; FLJ14466; Orai 1; ORAI calcium release activated calcium modulator 1; orai1; ORAT1; Protein orai 1; Protein orai-1; TMEM 142A; TMEM142A; Transmembrane protein 142A.</p>	<p>Isotype: IgG</p> <p>SWISS: Q96D31</p>	<p>Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ICC/IF (1:100-500) ELISA (1:5000-10000)</p> <p>Reactivity: (predicted: Human, Mouse, Rat, Rabbit, Pig, Sheep, Cow, Dog, Horse)</p> <p>Predicted MW.: 33 kDa</p> <p>Subcellular Location: Cell membrane</p>
--	--	---