
DFFA Antibody Blocking Peptide

Catalog Number: bs-0214P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Caspase-activated deoxyribonuclease(CAD) also designated DNA fragmentation factor 45 kDa subunit (DFF-45)and Caspase-activated DNase. CAD is a nuclease that induces DNA fragmentation and chromatin condensation during apoptosis. Degrades naked DNA and induces apoptotic morphology. ICAD (DNA fragmentation factor alpha subunit) (DNA fragmentation factor 45 kDa subunit) (DFF-45) (Inhibitor of CAD) is inhibitor of the caspase-activated DNase (DFF44).