

D-Dimer Mouse mAb

Cat. No: V1101
Name: D-Dimer Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: G1E9
Host/Source: Mouse
Subtype: IgG2b
Reactivity: (predicted: Human)
Format: Liquid
Concentration: $\geq 1\text{mg/ml}$
Buffer: 10mM PBS, pH7.4.
Applications: ELISA
Antigen: Native D-dimer protein derived from human plasma ([bs-41172P](#))
Pairing Information:

| ELISA | | Detection | |
|---------|-------|-----------|-------|
| | | V1101 | V1102 |
| Capture | V1101 | - | - |
| | V1102 | +++ | - |