

---

## **GPR4 Antibody Blocking Peptide**

Catalog Number: bs-13522P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background:

- Enables G protein-coupled receptor activity. Involved in several processes, including G protein-coupled receptor signaling pathway; positive regulation of Rho protein signal transduction; and response to acidic pH. Located in plasma membrane.  
[provided by Alliance of Genome Resources, Apr 2022]