
Recombinant CKC-1022 protein, Trx & S & His

Catalog Number: bs-49106P

Concentration: >0.5 mg/ml

AA Seq: 19-213/213

Predicted MW: 40

Tags: Trx & S & His

Activity: Not tested

Endotoxin: Not analyzed

Purity: >95% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris-HCl (pH8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.