
Recombinant human DC cocktail protein (Active, HEK293)

Catalog Number: bs-48032P

Concentration: >0.5 mg/ml

AA Seq: N/A

Tags: Tag free

Activity: Yes

Endotoxin: <0.5 EU/μg as determined by LAL.

Purity: >95% as determined by SDS-PAGE

Purification: AC

Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. In general: 12 months from date of receipt, -20 to -80°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 3 months, -20 to -80°C under sterile conditions after reconstitution.

Background: The cocktail contained recombinant human GM-CSF and IL-4. The concentrations of both cytokines had been optimized for dendritic cell culture.