
Recombinant human EPO protein (Active, HEK293)

Catalog Number: bs-48036P

Concentration: >0.5 mg/ml

AA Seq: 1-193/193

Detected MW: 54-67 kDa

Tags: Tag free

Activity: Yes

Endotoxin: <0.5 EU/μg as determined by LAL.

Purity: >95% as determined by SDS-PAGE

Purification: AC

Storage: Use a manual defrost freezer and avoid repeated freeze-thaw cycles. In general: 12 months from date of receipt, -20 to -80°C as supplied. 1 month, 2 to 8°C under sterile conditions after reconstitution. 3 months, -20 to -80°C under sterile conditions after reconstitution.

Background: Human erythropoietin is member of the EPO/TPO family and encodes a secreted, glycosylated cytokine hormone composed of four alpha helical bundles. The protein is found in the plasma and regulates red cell production by promoting erythroid differentiation and initiating hemoglobin synthesis. This protein also has neuroprotective activity against a variety of potential brain injuries and antiapoptotic functions in several tissue types. It is produced by kidney or liver of adult mammals and by liver of fetal or neonatal mammals.