

bsm-52676R**[Primary Antibody]**

www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

TFG Recombinant Rabbit mAb**— DATASHEET —**

Host: Rabbit	Isotype: IgG	Applications: WB (1:500-1000) IHC-P (1:100-500) IHC-F (1:400-800) IF (1:50-200) ICC/IF (1:50-200) Reactivity: (predicted: Human, Mouse, Rat) Predicted MW.: 43 kDa Subcellular Location: Nucleus
Clonality: Recombinant	CloneNo.: 5C1	
GeneID: 10342	SWISS: Q92734	
Target: TFG		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 1*TBS (pH7.4), 0.05% BSA, 40% Glycerol. Preservative: 0.02% Proclin300. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
Background: There are several documented fusion oncoproteins encoded partially by this gene. This gene also participates in several oncogenic rearrangements resulting in anaplastic lymphoma and mixoid chondrosarcoma, and may play a role in the NF-kappaB pathway. Multiple transcript variants have been found for this gene. [provided by RefSeq, Sep 2010]		