
TFG Recombinant Rabbit mAb

Catalog Number: bsm-52676R

Target Protein: TFG

Concentration: 1mg/ml

Form: Liquid

Host: Rabbit

Clonality: Recombinant

Clone No.: 5C1

Isotype: IgG

Applications: WB (1:500-1000), IHC-P (1:100-500), IHC-F (1:400-800), IF (1:50-200), ICC/IF (1:50-200)

Reactivity: (predicted:Human, Mouse, Rat)

Predicted MW: 43 kDa

Entrez Gene: 10342

Swiss Prot: Q92734

Purification: affinity purified by Protein A

Storage: 1*TBS (pH7.4), 0.05% BSA, 40% Glycerol. Preservative: 0.02% Proclin300.

Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: There are several documented fusion oncoproteins encoded partially by this gene. This gene also participates in several oncogenic rearrangements resulting in anaplastic lymphoma and mixoid chondrosarcoma, and may play a role in the NF-kappaB pathway. Multiple transcript variants have been found for this gene. [provided by RefSeq, Sep 2010]