

---

## **AKAP14 Antibody Blocking Peptide**

Catalog Number: bs-23424P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a member of A-kinase anchoring proteins (AKAPs), a family of functionally related proteins that target protein kinase A to discrete locations within the cell. The human ortholog of this gene has been characterized. It is reported to be enriched in airway cilia and to regulate ciliary beat frequency by anchoring the protein kinase A holoenzyme near substrates in the axoneme. [provided by RefSeq, May 2013].