

bs-4642R

[Primary Antibody]

KRV-VP1 + KRV-VP2 (C-terminus) Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

— DATASHEET —

Host: Rabbit	Isotype: IgG	Applications: IHC-P (1:100-500) IHC-F (1:100-500) IF (1:100-500) ELISA (1:5000-10000)
Clonality: Polyclonal		Reactivity: (predicted: Kilham virus)
Target: KRV-VP1 + KRV-VP2 (C-terminus)		Predicted MW.: 64 kDa
Immunogen: KLH conjugated synthetic peptide derived from Ki Kilham Rat Virus/KRV-VP1/KRV-VP2: 501-582/582.		
Purification: affinity purified by Protein A		
Concentration: 1mg/ml		
Storage: 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		

— SELECTED CITATIONS —

- **[IF=3.28]** Alkanani, Aimon K., et al. "Kilham rat virus-induced type 1 diabetes involves beta cell infection and intra-islet JAK-STAT activation prior to insulinitis." Virology 468 (2014): 19-27. IHC ;="Rat". 25129435