

Recombinant human ZO-1 protein, His

Catalog Number: bs-41327P

Concentration: >0.5mg/ml

AA Seq: 1467-1748/1748

Predicted MW: 33.2

Tags: N-His

Activity: Not tested

Purity: >90% as determined by SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 20mM Tris (pH8.0).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: This gene encodes a member of the membrane-associated guanylate kinase (MAGUK) family of proteins, and acts as a tight junction adaptor protein that also regulates adherens junctions. Tight junctions regulate the movement of ions and macromolecules between endothelial and epithelial cells. The multidomain structure of this scaffold protein, including a postsynaptic density 95/disc-large/zona occludens (PDZ) domain, a Src homology (SH3) domain, a guanylate kinase (GuK) domain and unique (U) motifs all help to co-ordinate binding of transmembrane proteins, cytosolic proteins, and F-actin, which are required for tight junction function. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Aug 2017]

VALIDATION IMAGES



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.