
Streptococcus suis 2 protein

Catalog Number: bs-4536P

AA Seq: Purified native protein

Activity: Not tested

Endotoxin: Not analyzed

Form: Lyophilized or Liquid

Storage: 10 mM TBS (pH=7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: Streptococcus suis serotype type 2 (SS2) infection has become a public health concern due to its zoonotic capability to cause severe infections in slaughterhouse workers and those who handle infected pork.