

---

## Recombined ASFV p72 protein

Catalog Number: bs-41384P

Concentration: >0.5 mg/ml

AA Seq: 1-646/646

Predicted MW: 73.1

Detected MW: 86 kDa

Tags: N-His Tag

Activity: No tested

Endotoxin: Not analyzed

Purity: >90% SDS-PAGE

Purification: AC

Form: Lyophilized or Liquid

Storage: 8M urea 1mM DTT, pH9.0.

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

**Background:** Major capsid protein, one of the most immunogenic protein in natural infections. Plays an essential role in virion assembly (PubMed:9580160).

### VALIDATION IMAGES

---



The purity of the protein is greater than 90% as determined by reducing SDS-PAGE.