
E.coli DH-5 Alpha Bacterial protein

Catalog Number: bs-2033P

AA Seq: Purified native protein

Activity: Not tested

Endotoxin: Not analyzed

Form: Lyophilized or Liquid

Storage: PBS (pH=7.4).

Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.

Background: This antibody reacts with Escherichia coli DH5a (E.coli DH5a), Escherichia coli is a gram negative bacillus that belongs to a larger group of Enterobacteriaceae - bacteria that inhabit the gastrointestinal tract. Although usually a harmless resident of the gut, some strains have the potential to cause serious problems, especially where there is an immature immune system or immunosuppression, or where the subtype of organism has acquired the ability to produce pathogenic toxins.

PRODUCT SPECIFIC PUBLICATIONS

[IF=4.916] Li X et al. A novel and facile immunosensor based on a barometer: Application for rapid analysis of Escherichia coli in waters. Talanta, 120859 (2020). Other ; . doi:10.1016/j.talanta.2020.120859

[IF=3.265] Jia X et al. Down-regulation of GCLC is involved in microcystin-LR-induced malignant transformation of human liver cells. Toxicology. 2019 Mar 30;421:49-58. Other ; . 30940545