

**bs-12700R**

**[ Primary Antibody ]**

## PLAC1 Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

|  |                      |   |
|--|----------------------|---|
| <b>Host:</b> Rabbit  | <b>Isotype:</b> IgG  | <b>Applications:</b> IHC-P (1:100-500)<br>IHC-F (1:100-500)<br>IF (1:100-500)<br>ICC/IF (1:100-500)<br>ELISA (1:5000-10000)<br><br><b>Reactivity:</b> (predicted: Human)<br><br><br><b>Predicted MW.:</b> 21 kDa<br><br><b>Subcellular Location:</b> Secreted |
| <b>Clonality:</b> Polyclonal   |                      |   |
| <b>GeneID:</b> 10761   | <b>SWISS:</b> Q9HBJ0 |   |
| <b>Target:</b> PLAC1   |                      |   |
| <b>Immunogen:</b> KLH conjugated synthetic peptide derived from human PLAC1: 101-180/212.  |                      |   |
| <b>Purification:</b> affinity purified by Protein A  |                      |   |
| <b>Concentration:</b> 1mg/ml   |                      |   |
| <b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.<br>Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles. |                      |   |