
SPATC1 Antibody Blocking Peptide

Catalog Number: bs-17636P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: SPATA15 is a 591 amino acid cytoplasmic protein expressed in testis that is part of a complex that includes p53 CDC, Cdc27. The gene encoding Speriolin is located on human chromosome 8, which consists of nearly 146 million base pairs, houses more than 800 genes and is associated with a variety of diseases and malignancies. Schizophrenia, bipolar disorder, Trisomy 8, Pfeiffer syndrome, congenital hypothyroidism, Waardenburg syndrome and some leukemias and lymphomas are thought to occur as a result of defects in specific genes that map to chromosome 8.