

---

## VIP Antibody Blocking Peptide

Catalog Number: bs-0077P

Activity: Yes

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: VIP(Vasoactive Intestinal peptide) is a 28-amino acid peptide structurally related to secretin.

It was originally isolated from intestinal extracts and shown to be a potent vasodilator.

Subsequent work demonstrated that VIP is very widely distributed in the peripheral and central nervous systems, and probably should not be considered a true GI hormone. A huge number of biological effects have been attributed to VIP. With respect to the digestive system, VIP seems to induce smooth muscle relaxation (lower esophageal sphincter, stomach, gallbladder), stimulate secretion of water into pancreatic juice and bile, and cause inhibition of gastric acid secretion and absorption from the intestinal lumen.