

## Human PRL Mouse mAb

Cat. No: V6602  
Name: Human PRL Mouse mAb  
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.  
Clone type: Monoclonal  
Clone No: 1D10  
Host/Source: Mouse  
Subtype: IgG1  
Reactivity: (predicted:&nbsp;Human)  
Format: Liquid  
Buffer: 0.01M PBS (pH7.4).  
Applications: ICA  
Antigen: Recombinant Human PRL Protein ([bs-43099P](#))  
Pairing Information:

| ELISA   |       | Detection |       |       |       |       |       |       |
|---------|-------|-----------|-------|-------|-------|-------|-------|-------|
|         |       | V6601     | V6602 | V6603 | V6604 | V6605 | V6606 | V6607 |
| Capture | V6601 | -         | -     | -     | -     | -     | -     | -     |
|         | V6602 | -         | -     | -     | -     | -     | -     | -     |
|         | V6603 | -         | -     | -     | -     | +++   | -     | -     |
|         | V6604 | -         | -     | -     | -     | -     | -     | -     |
|         | V6605 | -         | -     | -     | -     | -     | +++   | -     |
|         | V6606 | -         | -     | -     | -     | -     | -     | -     |
|         | V6607 | -         | -     | -     | -     | -     | +     | -     |
| ICA     |       | Detection |       |       |       |       |       |       |
|         |       | V6601     | V6602 | V6603 | V6604 | V6605 | V6606 | V6607 |
| Capture | V6601 | -         | -     | +++   | -     | -     | -     | -     |
|         | V6602 | -         | -     | +     | +++   | -     | -     | -     |
|         | V6603 | -         | -     | -     | +++   | -     | -     | -     |
|         | V6604 | -         | -     | -     | -     | -     | -     | -     |
|         | V6605 | -         | -     | -     | -     | -     | -     | -     |
|         | V6606 | -         | -     | -     | -     | -     | -     | -     |
|         | V6607 | -         | -     | -     | -     | -     | -     | -     |