

Human PRL Mouse mAb

Cat. No:	V6603
Name:	Human PRL Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	3F10
Host/Source:	Mouse
Subtype:	IgG1
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Recombinant Human PRL Protein (bs-43099P)
Pairing Information:	

ELISA		Detection						
		V6601	V6602	V6603	V6604	V6605	V6606	V6607
Capture	V6601	-	-	-	-	-	-	-
	V6602	-	-	-	-	-	-	-
	V6603	-	-	-	-	+++	-	-
	V6604	-	-	-	-	-	-	-
	V6605	-	-	-	-	-	+++	-
	V6606	-	-	-	-	-	-	-
	V6607	-	-	-	-	-	+	-
ICA		Detection						
		V6601	V6602	V6603	V6604	V6605	V6606	V6607
Capture	V6601	-	-	+++	-	-	-	-
	V6602	-	-	+	+++	-	-	-
	V6603	-	-	-	+++	-	-	-
	V6604	-	-	-	-	-	-	-
	V6605	-	-	-	-	-	-	-
	V6606	-	-	-	-	-	-	-
	V6607	-	-	-	-	-	-	-