

Human PRL Mouse mAb

Cat. No: V6606
Name: Human PRL Mouse mAb
Storage: Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type: Monoclonal
Clone No: 5F6
Host/Source: Mouse
Subtype: IgG1
Reactivity: (predicted: Human)
Format: Liquid
Buffer: 0.01M PBS (pH7.4).
Applications: ELISA
Antigen: Recombinant Human PRL Protein ([bs-43099P](#))
Pairing Information:

ELISA		Detection						
		V6601	V6602	V6603	V6604	V6605	V6606	V6607
Capture	V6601	-	-	-	-	-	-	-
	V6602	-	-	-	-	-	-	-
	V6603	-	-	-	-	+++	-	-
	V6604	-	-	-	-	-	-	-
	V6605	-	-	-	-	-	+++	-
	V6606	-	-	-	-	-	-	-
	V6607	-	-	-	-	-	+	-
ICA		Detection						
		V6601	V6602	V6603	V6604	V6605	V6606	V6607
Capture	V6601	-	-	+++	-	-	-	-
	V6602	-	-	+	+++	-	-	-
	V6603	-	-	-	+++	-	-	-
	V6604	-	-	-	-	-	-	-
	V6605	-	-	-	-	-	-	-
	V6606	-	-	-	-	-	-	-
	V6607	-	-	-	-	-	-	-