

Human Ferritin Light Chain Mouse mAb

Cat. No:	V2406
Name:	Human Ferritin Light Chain Mouse mAb
Storage:	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Clone type:	Monoclonal
Clone No:	8F7
Host/Source:	Mouse
Subtype:	IgG2b
Reactivity:	(predicted: Human)
Format:	Liquid
Buffer:	0.01M PBS (pH7.4).
Applications:	ELISA, ICA
Antigen:	Recombinant Human Ferritin light chain Protein (bs-41189P)
Pairing Information:	

ELISA		Detection			
		V2403	V2404	V2405	V2406
Capture	V2403	-	-	-	-
	V2404	-	-	+++	-
	V2405	-	+++	-	-
	V2406	-	-	+	-
ICA					
		V2403	V2404	V2405	V2406
Capture	V2403	-	+++	-	+++
	V2404	-	-	-	-
	V2405	-	-	-	-
	V2406	-	-	-	-