
CALCRL Antibody Blocking Peptide

Catalog Number: bs-1860P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: bs-1860P is one synthetic peptide derived from human CGRPR1.

CRLR is a receptor for calcitonin gene related peptide type 1. The activity of CRLR is mediated by G proteins which activate adenylyl cyclase. CRLR expression has been reported in brain, lung, blood vessel, liver, and intestinal tract. ESTs have been isolated from B Cell/lung/testis, bone marrow, embryo, lung, and synovium libraries.