

---

## **ZNF764 Antibody Blocking Peptide**

Catalog Number: bs-8661P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: ZNF764 (Zinc Finger Protein 764) is a Protein Coding gene. Diseases associated with ZNF764 include 16p11.2 microdeletion and thyroid hormone resistance. Among its related pathways are Gene Expression and Gene Expression. An important paralog of this gene is ZNF785.