

---

## PTPLAD1 Antibody Blocking Peptide

Catalog Number: bs-8693P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: PTPLAD1 is a component of Rac1-signaling pathways leading to the modulation of gene expression. It potentiates different effects of the small GTPase Rac1, such as c-Jun N-terminal kinase activation and transcriptional activity of nuclear factor kappaB (NF-kappaB). It plays a crucial role in HCV RNA replication and the propagation of JFH1 virus through interaction with viral and host proteins.