
SEH1L Antibody Blocking Peptide

Catalog Number: bs-8721P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: Seh1 is a 421 amino acid protein belonging to the WD repeat Sec13 family. Localized to the nucleus, Seh1 is a component of the nuclear pore complex Nup107-160. Nuclear pore complexes control bidirectional transport of macromolecules between the cytoplasm and the nucleus. All components of the complex Nup107-160, including Seh1, localize to the kinetochores during mitosis. Seh1 is expressed as two isoforms produced by alternative splicing and contains six WD repeats.