

---

## **SYNPO2 Antibody Blocking Peptide**

Catalog Number: bs-8743P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: SYNPO2 (Synaptopodin 2) is a Protein Coding gene. Diseases associated with SYNPO2 include Duchenne Muscular Dystrophy and Myopathy, Myofibrillar, 2. GO annotations related to this gene include actin binding and muscle alpha-actinin binding. An important paralog of this gene is SYNPO2L.