
KLF3 Antibody Blocking Peptide

Catalog Number: bs-8854P

Activity: Not tested

Purification: HPLC

Storage: Shipped at 4°C. Stored at -20°C for one year. Avoid repeated freeze/thaw cycles.

Background: KLF3 is a transcriptional repressor and activator. It also binds the GT-box of cyclin D1 promoter and mediates cell cycle progression at G phase as a downstream target of focal adhesion kinase (FAK). Ubiquitously expressed. Kruppel-like factor proteins all contain a characteristic C-terminal DNA binding domain containing 3 Kruppel-like zinc fingers, and they all recognize CACCC and GC boxes present in promoters and enhancers.