

**bs-17759R**

**[ Primary Antibody ]**

## MRGPRH Rabbit pAb



www.bioss.com.cn

sales@bioss.com.cn

techsupport@bioss.com.cn

400-901-9800

### — DATASHEET —

<b>Host:</b> Rabbit	<b>Isotype:</b> IgG	<b>Applications:</b> WB (1:500-2000)
<b>Clonality:</b> Polyclonal		<b>Reactivity:</b> (predicted: Rat)
<b>GeneID:</b> 404641	<b>SWISS:</b> Q7TN38	
<b>Target:</b> MRGPRH		<b>Predicted MW.:</b> 37 kDa
<b>Immunogen:</b> KLH conjugated synthetic peptide derived from rat MRGPRH: 1-100/321. < Extracellular >		<b>Subcellular Location:</b> Cell membrane
<b>Purification:</b> affinity purified by Protein A		
<b>Concentration:</b> 1mg/ml		
<b>Storage:</b> 0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol. Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.		
<b>Background:</b> seven transmembrane domain G-protein coupled receptor; most often found in sensory neurons [RGD, Feb 2006]		